The Board of Pharmacy has received notice of the following product recall. The Board strongly encourages pharmacies to immediately review their quality assurance and recall policies and procedures to determine if any corrective action is required.

Product	SKU	Lot
Gown Isolation W/Ties Univ Ylw	1100PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties Univ Ylw	V1100PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties XI Ylw	1101PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties Univ Blu	1200PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties Univ Ylw	2100PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties Univ Ylw	V2100PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties Univ Ylw	2150PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties Univ Blue	2200PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties XI Blue	2201PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Lab Coats SMS Large Blue	2202LC	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Lab Coats SMS X-Large Blue	2203LC	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Lab Coats SMS 2X-Large Blue	2204LC	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Lab Coats SMS Medium White	2301LC	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Lab Coats SMS Large White	2302LC	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Lab Coats SMS 2X-Large White	2304LC	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties Univ Ylw	V2100PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Isolation W/Ties Univ Blue	552200PG	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX

Product	SKU	Lot
Gown Isolation W/Ties Univ Ylw	N()N/4/58N	All lot numbers from <b>18K</b> X <b>B</b> XXXXX to <b>20K</b> X <b>B</b> XXXXX
Gown Kendall Blue RG	CT5000	All lot numbers from <b>1901</b> XX <b>K</b> to <b>2010</b> XX <b>K</b>
Gown Kendall Blue L	CT5100	All lot numbers from <b>1901</b> XX <b>K</b> to <b>2010</b> XX <b>K</b>
Gown Kendall Blue L	CT5101	All lot numbers from <b>1901</b> XX <b>K</b> to <b>2010</b> XX <b>K</b>
Gown Kendall Blue XXL	CT5102	All lot numbers from <b>1901</b> XX <b>K</b> to <b>2010</b> XX <b>K</b>

Cardinal Health has issued a Medical Device Recall for specific production lots of **Cardinal Health™ Isolation Gowns, Lab Coats and Kendall™ Gowns** due to potential rips, holes, tears, incomplete and/or open seams, staining, foreign material or debris on the gowns and lab coats.